



K-REX™ RONGEUR

SPECIALIST



K-Rex™ Rongeur



GEISTER MEDIZINTECHNIK GMBH

Föhrenstrasse 2 | D-78532 Tuttlingen | Germany

☎ +49 7461/966240 | ☎ +49 7461/9662422

www.geister.com | info@geister.com

“

The beginning of an end is sometimes just the start of something new.

And once in a while, it's the genesis of something wonderful.

— Jaz Primo



Fingertip design curved convex grip for optimum control.



A pleasure to HANDLE and TOUCH ...


GEISTER®

Fairgrip™ultra Handle Design

The Fairgrip™ultra ergonomic instrument handles with a new dimple knurling similar to the design of golfballs offer a completely new experience of precision, sensitivity and control during operations and help to reduce fatigue.



Our new handle design follows the convex hull of hand contour for a more comfortable grip, less fatigue and full control as every finger can grab the handle safely. This new and unique handle design is pleasure to touch and to work with.



Concave

Convex

A pleasure to HANDLE and TOUCH ...

7
GEISTER®





patent pending



A pleasure to MAINTAIN ...




Clickable conversion
open & close

- 👍 No parts can be lost
- 👍 No extra check for components
- 👍 No tools needed
- 👍 So easy!!!

Instruction for cleaning:



Ready for use:



A pleasure to MAINTAIN ...

GEISTER®

Convertible Style

By pushing the front handle forward the Kerrison punch is just ready for cleaning and lubrication - smart decoupling without losing any parts.






A pleasure to USE ...



Perfectly matched



Strong fixation points

A pleasure to USE ...

GEISTER®





K-Rex™ Rongeur



For more
information, visit:



A pleasure to HANDLE and TOUCH ...

Experience exceptional comfort, safety, and control with the ergonomically designed Fairgrip™ultra handle. Its refined texture and one-size-fits-all design ensure a secure, enjoyable grip and a pleasant tactile experience for every user.

to MAINTAIN ...

Simplify your workflow with the tool's quick and easy decoupling mechanism, designed for hassle-free cleaning and maintenance without the risk of losing parts.

to USE ...

Delivering unmatched performance, the precisely engineered components provide a powerful and precise bite, ensuring reliability and innovation in every procedure.




K-Rex™ Rongeur

For a precise removal of bone and tissue in General Surgery,
Traumatology, Orthopaedics.



40-3220C3




Golden handle indicates
thin footplate

40-3320C3



SurgiBOX™ Micro

duplex with 2 permanent filters
for max. 8 K-Rex™ punches

Art.-No.	Measure mm
77-0120.01	495 × 270 × 75





Rongeurs not included



Example

77-4150



Example

77-4151


Storage Tray Basket for K-Rex

for Surgibox™ 450, 3/4 size


Art.-No.	Measure mm
77-4150	396 x 250 x 65

for Surgibox™ 600, 1/1 size

Art.-No.	Measure mm
774151	539 x 250 x 65



Regular footplate




40° up


K-Rex™ Rongeur

detachable

Length	Art.-No. Width 1 mm without ejector	Art.-No. Width 2 mm with ejector	Art.-No. Width 3 mm with ejector	Art.-No. Width 4 mm with ejector	Art.-No. Width 5 mm with ejector	Art.-No. Width 6 mm with ejector
180 mm	40-3218C1	40-3218C2	40-3218C3	40-3218C4	40-3218C5	40-3218C6
200 mm	40-3220C1	40-3220C2	40-3220C3	40-3220C4	40-3220C5	40-3220C6
230 mm	40-3223C1	40-3223C2	40-3223C3	40-3223C4	40-3223C5	40-3223C6
280 mm	40-3228C1	40-3228C2	40-3228C3	40-3228C4	40-3228C5	40-3228C6



Thin footplate*



40° up

180 mm	40-3318C1	40-3318C2	40-3318C3	40-3318C4	40-3318C5	40-3318C6
200 mm	40-3320C1	40-3320C2	40-3320C3	40-3320C4	40-3320C5	40-3320C6
230 mm	40-3323C1	40-3323C2	40-3323C3	40-3323C4	40-3323C5	40-3323C6

The better way to operate™

GEISTER®

78532 Tuttlingen, Germany | www.geister.com

